

Antibodies Product Datasheet

Product Name:	CNKS1 Antibody
Catalog Number:	R12-2634
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	CNKS1 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CNKS1.
Concentration:	1 mg/ml
Gene Name:	CNKS1
Other Name:	Connector enhancer of kinase suppressor of ras 1; Connector enhancer of KSR 1; CNK homolog protein 1; CNK1; hCNK1; Connector enhancer of KSR-like; CNKS1; CNK1
Human Swiss Prot No.:	Q969H4
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	CNKS1 Antibodyantibody detects endogenous levels of CNKS1.
Storage and Stability:	Stable at -20°C for at least 1 year.