

Antibodies Product Datasheet

Product Name:	IL-2R α Antibody
Catalog Number:	C30337
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human IL2RA. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	IL2RA
Protein Name:	Interleukin-2 receptor subunit alpha
Other Name:	IL2RA; Interleukin-2 receptor subunit alpha; IL-2 receptor subunit alpha; IL-2-RA; IL-2R subunit alpha; IL2-RA; TAC antigen; p55; CD25
Human Swiss Prot No.:	P01589
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-2R α Polyclonal Antibody detects endogenous levels of IL-2R α protein.
Storage and Stability:	-20°C/1 year