

## Antibodies Product Datasheet

<b>Product Name:</b>	c-Fms (Phospho-Tyr708) Antibody
<b>Catalog Number:</b>	A9015
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized phospho-peptide around the phosphorylation site of human c-Fms (phospho Tyr708)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	CSF1R
<b>Protein Name:</b>	Macrophage colony-stimulating factor 1 receptor
<b>Other Name:</b>	CSF1R; FMS; Macrophage colony-stimulating factor 1 receptor; CSF-1 receptor; CSF-1-R; CSF-1R; M-CSF-R; Proto-oncogene c-Fms; CD antigen CD115
<b>Human Swiss Prot No.:</b>	P07333
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Phospho-c-Fms (Y708) Polyclonal Antibody detects endogenous levels of c-Fms protein only when phosphorylated at Y708.
<b>Storage and Stability:</b>	-20°C/1 year