

---

## Antibodies Product Datasheet

<b>Product Name:</b>	p53 Monoclonal Antibody(6C4)
<b>Catalog Number:</b>	N1317
<b>Reactivity:</b>	Human
<b>Applications:</b>	WB,IHC-p,IF/ICC
<b>Dilutions:</b>	WB: 1:2000 IHC: 1:200 IF: 1:100-200
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthetic Peptide of p53
<b>Concentration:</b>	
<b>Gene Name:</b>	TP53
<b>Protein Name:</b>	Cellular tumor antigen p53
<b>Other Name:</b>	TP53; P53; Cellular tumor antigen p53; Antigen NY-CO-13; Phosphoprotein p53; Tumor suppressor p53
<b>Human Swiss Prot No.:</b>	P04637
<b>Formulation:</b>	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Specificity:</b>	The antibody detects endogenous p53 proteins.
<b>Storage and Stability:</b>	-20°C/1 year