

Antibodies Product Datasheet

Product Name:	ErbB-3 Monoclonal Antibody
Catalog Number:	N0240
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of ErbB-3 (aa1175-1275) expressed in E. Coli.
Concentration:	
Gene Name:	ERBB3
Protein Name:	Receptor tyrosine-protein kinase erbB-3
Other Name:	ERBB3; HER3; Receptor tyrosine-protein kinase erbB-3; Proto-oncogene-like protein c-ErbB-3; Tyrosine kinase-type cell surface receptor HER3
Human Swiss Prot No.:	P21860
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	ErbB-3 Monoclonal Antibody detects endogenous levels of ErbB-3 protein.
Storage and Stability:	-20°C/1 year