

Antibodies Product Datasheet

Product Name:	DHDDS Antibody
Catalog Number:	C50240
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DHDDS
Concentration:	1 mg/ml
Gene Name:	DHDDS HDS
Protein Name:	Dehydrolipoyl diphosphate synthase (Dedol-PP synthase) (EC 2.5.1.-) (Cis-isoprenyltransferase) (CIT) (Cis-IPTase) (Epididymis tissue protein Li 189m)
Other Name:	
Human Swiss Prot No.:	Q86SQ9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DHDDS at Human, Mouse
Storage and Stability:	Store at -20° C/1 year