

Antibodies Product Datasheet

Product Name:	Actin α 3 Antibody
Catalog Number:	C0123
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Actin-gamma2. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	ACTG2
Protein Name:	Actin gamma-enteric smooth muscle
Other Name:	ACTG2; ACTA3; ACTL3; ACTSG; Actin; gamma-enteric smooth muscle; Alpha-actin-3; Gamma-2-actin; Smooth muscle gamma-actin
Human Swiss Prot No.:	P63267
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Actin α 3 Polyclonal Antibody detects endogenous levels of Actin α 3 protein.
Storage and Stability:	-20°C/1 year