
Antibodies Product Datasheet

Product Name:	GRIK4 Antibody
Catalog Number:	C44003
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GRIK4
Concentration:	1 mg/ml
Gene Name:	GRIK4 GRIK
Protein Name:	GRIK4
Other Name:	Glutamate receptor, ionotropic kainate 4 (Excitatory amino acid receptor 1) (EAA1) (Glutamate receptor KA-1) (KA1)
Human Swiss Prot No.:	Q16099
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.239% sodium azide.
Specificity:	This antibody detects endogenous levels of GRIK4 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year