

Antibodies Product Datasheet

Product Name:	MEK1 Antibody
Catalog Number:	R12-3013
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	MEK1 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MEK1.
Concentration:	1 mg/ml
Gene Name:	MAP2K1
Other Name:	Dual specificity mitogen-activated protein kinase kinase 1; ERK activator kinase 1; MAP kinase kinase 1; MAP2K1; MAPK/ERK kinase 1; MAPKK 1; MAPKK1; MEK1; MP2K1; PRKMK1; kinase MEK1
Human Swiss Prot No.:	Q02750
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	MEK1 Antibody antibody detects endogenous levels of MEK1.
Storage and Stability:	Stable at -20°C for at least 1 year.