

Antibodies Product Datasheet

Product Name:	Serine/threonine protein kinase UL97 Antibody
Catalog Number:	C41530
Reactivity:	Human cytomegalovirus (strain AD169) (HHV-5)
Applications:	ELISA
Dilutions:	ELISA(peptide) 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Serine/threonine protein kinase UL97
Concentration:	1 mg/ml
Gene Name:	UL97
Protein Name:	Serine/threonine protein kinase UL97
Other Name:	
Human Swiss Prot No.:	P16788
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Serine/threonine protein kinase UL97
Storage and Stability:	-20°C/1 year