
Antibodies Product Datasheet

Product Name:	Fatty Acid Synthase Antibody
Catalog Number:	R12-2140
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p: 100-300. Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Fatty Acid Synthase. AA range:1478-1527
Concentration:	1 mg/ml
Gene Name:	FASN
Protein Name:	Fatty acid synthase
Other Name:	FASN; FAS; Fatty acid synthase
Human Swiss Prot No.:	P49327
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Fatty Acid Synthase Polyclonal Antibody detects endogenous levels of Fatty Acid Synthase protein.
Storage and Stability:	-20°C/1 year