

Antibodies Product Datasheet

Product Name:	Anti-Human Cytokeratin 14 (ABT049) Mouse Monoclonal Antibody
Catalog Number:	V0290
Clone ID:	ABT049
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	KRT14
Gene Name:	P02533
Other Name:	CK 14, CK-14, ck14, Cytokeratin 14, Cytokeratin-14, Cytokeratin14, Dowling Meara, EBS3, EBS4, Epidermolysis bullosa simplex, K14, K1C14 HUMAN, Keratin 14 (epidermolysis bullosa simplex, Dowling-Meara, Koebner), Keratin 14, Keratin, Keratin type I cytoskeletal 14, Keratin, type I cytoskeletal 14, Keratin-14, Keratin14, Koebner, Krt 14, Krt14, NFJ, OTTHUMP00000164624, type I cytoskeletal 14
Isotype	IgG1, Kappa
Localization:	Tonsil/ Prostate
Positive Tissue Control:	Cytoplasmic, Membranous
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human CK14 protein. In immunohistochemistry on formalin-fixed, paraffin-embedded tissue sections, the antibody specifically labels the basal cell of squamous epithelial cells and glandular epithelia, myoepithelium a
Storage and Stability:	-20°C/1 year