

Antibodies Product Datasheet

Product Name:	NUDT3/4/10/11 Antibody
Catalog Number:	C49972
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human NUDT3/4/10/11
Concentration:	1 mg/ml
Gene Name:	NUDT3 DIPP DIPP1
Protein Name:	Diphosphoinositol polyphosphate phosphohydrolase 1 (DIPP-1) (EC 3.6.1.52) (Diadenosine 5',5'''-P1,P6-hexaphosphate hydrolase 1) (EC 3.6.1.-) (Nucleoside diphosphate-linked moiety X motif 3) (Nudix motif 3)
Other Name:	
Human Swiss Prot No.:	O95989
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of NUDT3/4/10/11 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year