

Antibodies Product Datasheet

Product Name:	Flt3-L Antibody
Catalog Number:	C30133
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human FLT3LG. AA range:171-220
Concentration:	1 mg/ml
Gene Name:	FLT3LG
Protein Name:	Fms-related tyrosine kinase 3 ligand
Other Name:	FLT3LG; Fms-related tyrosine kinase 3 ligand; Flt3 ligand; Flt3L; SL cytokine
Human Swiss Prot No.:	P49771
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Flt3-L Polyclonal Antibody detects endogenous levels of Flt3-L protein.
Storage and Stability:	-20°C/1 year