

---

## Antibodies Product Datasheet

<b>Product Name:</b>	COX1 Antibody
<b>Catalog Number:</b>	C47461
<b>Reactivity:</b>	Human
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB: 1:1000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthetic Peptide of COX1
<b>Concentration:</b>	
<b>Gene Name:</b>	PTGS1
<b>Protein Name:</b>	Prostaglandin G/H synthase 1
<b>Other Name:</b>	PTGS1; COX1; Prostaglandin G/H synthase 1; Cyclooxygenase-1; COX-1; Prostaglandin H2 synthase 1; PGH synthase 1; PGHS-1; PHS 1; Prostaglandin-endoperoxide synthase 1; YM3223
<b>Human Swiss Prot No.:</b>	P23219
<b>Formulation:</b>	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Specificity:</b>	The antibody detects endogenous COX1 proteins.
<b>Storage and Stability:</b>	-20°C/1 year