

Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | RFX1 Antibody |
| Catalog Number: | C51162 |
| Reactivity: | Human, Mouse |
| Applications: | IHC, WB |
| Dilutions: | WB 1:500-2000 IHC 1:50-200 |
| Modification: | Unmodified/Total |
| Source: | Polyclonal, Rabbit,IgG |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | Synthesized peptide derived from human N-terminal RFX1 |
| Concentration: | 1 mg/ml |
| Gene Name: | RFX1 |
| Protein Name: | MHC class II regulatory factor RFX1 (Enhancer factor C) (EF-C) (Regulatory factor X 1) (RFX) (Transcription factor RFX1) |
| Other Name: | MHC class II regulatory factor RFX1 (Enhancer factor C) (EF-C) (Regulatory factor X 1) (RFX) (Transcription factor RFX1) |
| Human Swiss Prot No.: | P22670 |
| Formulation: | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide. |
| Specificity: | This antibody detects endogenous levels of RFX1 at Human, Mouse |
| Storage and Stability: | Store at -20° C/1 year |