

Antibodies Product Datasheet

Product Name:	PIAS 1 Antibody
Catalog Number:	C0361
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PIAS1. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	PIAS1
Protein Name:	E3 SUMO-protein ligase PIAS1
Other Name:	PIAS1; DDXBP1; E3 SUMO-protein ligase PIAS1; DEAD/H box-binding protein 1; Gu-binding protein; GBP; Protein inhibitor of activated STAT protein 1; RNA helicase II-binding protein
Human Swiss Prot No.:	O75925
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PIAS 1 Polyclonal Antibody detects endogenous levels of PIAS 1 protein.
Storage and Stability:	-20°C/1 year