

Antibodies Product Datasheet

Product Name:	Cyclin D3 (Phospho-Thr283) Antibody
Catalog Number:	A0418
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Cyclin D3 around the phosphorylation site of Thr283. AA range:243-292
Concentration:	1 mg/ml
Gene Name:	CCND3
Protein Name:	G1/S-specific cyclin-D3
Other Name:	CCND3; G1/S-specific cyclin-D3
Human Swiss Prot No.:	P30281
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Cyclin D3 (T283) Polyclonal Antibody detects endogenous levels of Cyclin D3 protein only when phosphorylated at T283.
Storage and Stability:	-20°C/1 year