
Antibodies Product Datasheet

Product Name:	TPPP (Phospho-Ser18) Antibody
Catalog Number:	A8441
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human TPPP (Phospho-Ser18)
Concentration:	1 mg/ml
Gene Name:	TPPP TPPP1
Protein Name:	TPPP (Phospho-Ser18)
Other Name:	Tubulin polymerization-promoting protein (TPPP) (25 kDa brain-specific protein) (TPPP/p25) (p24) (p25-alpha)
Human Swiss Prot No.:	O94811
Formulation:	Liquid in PBS containing 50% glycerol, and 0.179% sodium azide.
Specificity:	This antibody detects endogenous levels of TPPP (Phospho-Ser18) at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year