

Antibodies Product Datasheet

Product Name:	Anti-Claudin 4 [ABT-CLD4]
Catalog Number:	V0061
Clone ID:	ABT-CLD4
Reactivity:	Human
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	O14493
Gene Name:	CLDN4
Other Name:	Claudin-4; CLD4; CLDN4; Clostridium perfringens enterotoxin receptor; CPE-receptor; CPE-R; CPER; CPETR1; WBSCR8; Williams-Beuren syndrome chromosome region 8 protein
Isotype	Mouse IgG
Localization:	N/A
Positive Tissue Control:	N/A
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	Anti-Claudin 4 [ABT-CLD4] detects endogenous levels of total Claudin 4 protein.
Storage and Stability:	-20°C/1 year