

Antibodies Product Datasheet

Product Name:	RBM22 Antibody
Catalog Number:	C50636
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human RBM22
Concentration:	1 mg/ml
Gene Name:	RBM22 ZC3H16 199G4
Protein Name:	Pre-mRNA-splicing factor RBM22 (RNA-binding motif protein 22) (Zinc finger CCCH domain-containing protein 16)
Other Name:	
Human Swiss Prot No.:	Q9NW64
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of RBM22 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year