

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | SR-2A Antibody |
| Catalog Number: | C12013 |
| Reactivity: | Human, Rat, Mouse |
| Applications: | IHC-p, IF/ICC, ELISA |
| Dilutions: | Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications. |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human 5-HT-2A. AA range:422-471 |
| Concentration: | 1 mg/ml |
| Gene Name: | HTR2A |
| Protein Name: | 5-hydroxytryptamine receptor 2A |
| Other Name: | HTR2A; HTR2; 5-hydroxytryptamine receptor 2A; 5-HT-2; 5-HT-2A; Serotonin receptor 2A |
| Human Swiss Prot No.: | P28223 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | SR-2A Polyclonal Antibody detects endogenous levels of SR-2A protein. |
| Storage and Stability: | -20°C/1 year |