

Antibodies Product Datasheet

Product Name:	MASP1 (light chain, Cleaved-Ile449) Antibody
Catalog Number:	L0408
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MASP1 (light chain, Cleaved-Ile449)
Concentration:	1 mg/ml
Gene Name:	MASP1 CRARF CRARF1 PRSS5
Protein Name:	MASP1 (light chain, Cleaved-Ile449)
Other Name:	Mannan-binding lectin serine protease 1 (EC 3.4.21.-) (Complement factor MASP-3) (Complement-activating component of Ra-reactive factor) (Mannose-binding lectin-associated serine protease 1) (MASP-1) (Mannose-binding protein-associated serine protease) (R
Human Swiss Prot No.:	P48740
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human MASP1 (light chain, Cleaved-Ile449, protein was cleaved amino acid sequence between 448-449)
Storage and Stability:	-20°C/1 year