

Antibodies Product Datasheet

Product Name:	MRP-L35 Antibody
Catalog Number:	C14073
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from MRP-L35 . at AA range: 80-160
Concentration:	1 mg/ml
Gene Name:	MRPL35
Protein Name:	39S ribosomal protein L35 mitochondrial
Other Name:	MRPL35; BM-007; 39S ribosomal protein L35, mitochondrial; L35mt; MRP-L35
Human Swiss Prot No.:	Q9NZE8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MRP-L35 Polyclonal Antibody detects endogenous levels of MRP-L35 protein.
Storage and Stability:	-20°C/1 year