

Antibodies Product Datasheet

Product Name:	NU3M Antibody
Catalog Number:	C45098
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 20-100
Concentration:	1 mg/ml
Gene Name:	MT-ND3 MTND3 NADH3 ND3
Protein Name:	NADH-ubiquinone oxidoreductase chain 3 (EC 1.6.5.3) (NADH dehydrogenase subunit 3)
Other Name:	
Human Swiss Prot No.:	P03897
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	NU3M Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year