

Antibodies Product Datasheet

Product Name:	Fes Antibody
Catalog Number:	C10294
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FES. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	FES
Protein Name:	Tyrosine-protein kinase Fes/Fps
Other Name:	FES; FPS; Tyrosine-protein kinase Fes/Fps; Feline sarcoma/Fujinami avian sarcoma oncogene homolog; Proto-oncogene c-Fes; Proto-oncogene c-Fps; p93c-fes
Human Swiss Prot No.:	P07332
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Fes Polyclonal Antibody detects endogenous levels of Fes protein.
Storage and Stability:	-20°C/1 year