

## Antibodies Product Datasheet

<b>Product Name:</b>	SLC6A15 Antibody
<b>Catalog Number:</b>	C17577
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human SLC6A15. AA range:32-81
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	SLC6A15
<b>Protein Name:</b>	Sodium-dependent neutral amino acid transporter B(0)AT2
<b>Other Name:</b>	SLC6A15; B0AT2; NTT73; SBAT1; Sodium-dependent neutral amino acid transporter B(0)AT2; Sodium- and chloride-dependent neurotransmitter transporter NTT73; Sodium-coupled branched-chain amino-acid transporter 1; Solute carrier family 6 member
<b>Human Swiss Prot No.:</b>	Q9H2J7
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	SLC6A15 Polyclonal Antibody detects endogenous levels of SLC6A15 protein.
<b>Storage and Stability:</b>	-20°C/1 year