

Antibodies Product Datasheet

Product Name:	ERI3 Antibody
Catalog Number:	C43639
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ERI3
Concentration:	1 mg/ml
Gene Name:	ERI3 PINT1 PRNPIP PRNPIP1
Protein Name:	ERI1 exoribonuclease 3 (EC 3.1.-.-) (Prion interactor 1) (Prion protein-interacting protein)
Other Name:	
Human Swiss Prot No.:	O43414
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ERI3
Storage and Stability:	-20°C/1 year