

Antibodies Product Datasheet

Product Name:	Glial Fibrillary Acidic Protein (GFAP) (PT1995) mouse mAb
Catalog Number:	N1844
Reactivity:	Human, Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC-p 1:100-500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse/IgG1, Kappa
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Glial Fibrillary Acidic Protein (GFAP)
Concentration:	
Gene Name:	GFAP
Protein Name:	Glial fibrillary acidic protein (GFAP)
Other Name:	
Human Swiss Prot No.:	P14136
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Glial Fibrillary Acidic Protein (GFAP). Heat-induced epitope retrieval (HIER) TRIS-EDTA of pH8.0 was highly recommended as antigen repair method in paraffin section
Storage and Stability:	-20°C/1 year