

Antibodies Product Datasheet

Product Name:	PCSK7 Antibody
Catalog Number:	C50108
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human PCSK7
Concentration:	1 mg/ml
Gene Name:	PCSK7 LPC PC7 PC8 SPC7
Protein Name:	Proprotein convertase subtilisin/kexin type 7 (EC 3.4.21.-) (Lymphoma proprotein convertase) (Prohormone convertase 7) (Proprotein convertase 7) (PC7) (Proprotein convertase 8) (PC8) (hPC8) (Subtilisin/kexin-like protease PC7)
Other Name:	
Human Swiss Prot No.:	Q16549
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PCSK7 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year