

Antibodies Product Datasheet

Product Name:	OR6C68 Antibody
Catalog Number:	R12-3250
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	OR6C68 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human OR6C68.
Concentration:	1 mg/ml
Gene Name:	OR6C68
Other Name:	Olfactory receptor 6C68; OR6C68
Human Swiss Prot No.:	A6NDL8
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	OR6C68 Antibodyantibody detects endogenous levels of OR6C68.
Storage and Stability:	Stable at -20°C for at least 1 year.