

Antibodies Product Datasheet

Product Name:	GAS2L1 Antibody
Catalog Number:	C51021
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GAS2L1
Concentration:	1 mg/ml
Gene Name:	GAS2L1 GAR22
Protein Name:	GAS2-like protein 1 (GAS2-related protein on chromosome 22) (Growth arrest-specific protein 2-like 1)
Other Name:	
Human Swiss Prot No.:	Q99501
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GAS2L1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year