

## Antibodies Product Datasheet

<b>Product Name:</b>	IPP-2 (Phospho-Ser44) Antibody
<b>Catalog Number:</b>	A8402
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized phospho-peptide around the phosphorylation site of human IPP-2 (phospho Ser44)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	PPP1R2
<b>Protein Name:</b>	Protein phosphatase inhibitor 2
<b>Other Name:</b>	PPP1R2; IPP2; Protein phosphatase inhibitor 2; IPP-2
<b>Human Swiss Prot No.:</b>	P41236
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Phospho-IPP-2 (S44) Polyclonal Antibody detects endogenous levels of IPP-2 protein only when phosphorylated at S44.
<b>Storage and Stability:</b>	-20°C/1 year