
Antibodies Product Datasheet

Product Name:	Cleaved-COL4A3 (L1425) Antibody
Catalog Number:	L0227
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Collagen IV alpha3. AA range:1376-1425
Concentration:	1 mg/ml
Gene Name:	COL4A3
Protein Name:	Collagen alpha-3(IV) chain
Other Name:	COL4A3; Collagen alpha-3(IV) chain; Goodpasture antigen
Human Swiss Prot No.:	Q01955
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cleaved-COL4A3 (L1425) Polyclonal Antibody detects endogenous levels of fragment of activated COL4A3 protein resulting from cleavage adjacent to L1425.
Storage and Stability:	-20°C/1 year