

Antibodies Product Datasheet

Product Name:	PDGFR α mouse Monoclonal Antibody(7A3)
Catalog Number:	N1662
Reactivity:	Human,Rat,Mouse
Applications:	IF/ICC,IHC-p
Dilutions:	IF: 1:50-200 IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of PDGFR α at AA range of 1010-1090
Concentration:	1 mg/ml
Gene Name:	PDGFRA
Protein Name:	PDGFRA
Other Name:	PDGFRA
Human Swiss Prot No.:	P16234
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PDGFR α protein detects endogenous levels of PDGFRA
Storage and Stability:	-20°C/1 year