

Antibodies Product Datasheet

Product Name:	COP ζ 1 Antibody
Catalog Number:	C15184
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human COPZ1. AA range:11-60
Concentration:	1 mg/ml
Gene Name:	COPZ1
Protein Name:	Coatomer subunit zeta-1
Other Name:	COPZ1; COPZ; CGI-120; HSPC181; Coatomer subunit zeta-1; Zeta-1-coat protein; Zeta-1 COP
Human Swiss Prot No.:	P61923
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	COP ζ 1 Polyclonal Antibody detects endogenous levels of COP ζ 1 protein.
Storage and Stability:	-20°C/1 year