

Antibodies Product Datasheet

Product Name:	WIP Antibody
Catalog Number:	C19552
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human WIPF1. AA range:421-470
Concentration:	1 mg/ml
Gene Name:	WIPF1
Protein Name:	WAS/WASL-interacting protein family member 1
Other Name:	WIPF1; WASPIP; WIP; WAS/WASL-interacting protein family member 1; Protein PRPL-2; Wiskott-Aldrich syndrome protein-interacting protein; WASP-interacting protein
Human Swiss Prot No.:	O43516
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	WIP Polyclonal Antibody detects endogenous levels of WIP protein.
Storage and Stability:	-20°C/1 year