

Antibodies Product Datasheet

Product Name:	TMP21 Antibody
Catalog Number:	C48413
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human TMP21
Concentration:	1 mg/ml
Gene Name:	TMED10 TMP21
Protein Name:	Transmembrane emp24 domain-containing protein 10 (21 kDa transmembrane-trafficking protein) (S31III125) (S31I125) (Tmp-21-I) (Transmembrane protein Tmp21) (p23) (p24 family protein delta-1) (p24delta1) (p24delta)
Other Name:	
Human Swiss Prot No.:	P49755
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of TMP21 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year