

Antibodies Product Datasheet

Product Name:	MAPKAPK5 (Phospho-Thr182) Antibody
Catalog Number:	P12-1144
Reactivity:	Human, Mouse
Applications:	WB, IHC, ELISA
Dilutions:	MAPKAPK5 (Phospho-Thr182) Antibody
Modification:	Phospho Specific
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using phospho peptide. The antibody against non-phospho peptide was removed by chromatography using corresponding non-phospho peptide.
Immunogen:	The antiserum was produced against synthesized phospho-peptide derived from human MAPKAPK5 around the phosphorylation site of Thr182.
Concentration:	1 mg/ml
Gene Name:	MAPKAPK5
Other Name:	MAP kinase-activated protein kinase 5; MAPK-activated protein kinase 5; MAPK5; MAPKAP kinase 5; p38-regulated/activated protein kinase; PRAK
Human Swiss Prot No.:	Q8IW41
Formulation:	Rabbit IgG in phosphate buffered saline (Phospho-without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	MAPKAPK5 (Phospho-Thr182) Antibody detects endogenous levels of MAPKAPK5 only when phosphorylated at Thr182.
Storage and Stability:	Stable at -20°C for at least 1 year.