

Antibodies Product Datasheet

Product Name:	NME6 Antibody
Catalog Number:	C50146
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human NME6
Concentration:	1 mg/ml
Gene Name:	NME6
Protein Name:	Nucleoside diphosphate kinase 6 (NDK 6) (NDP kinase 6) (EC 2.7.4.6) (Inhibitor of p53-induced apoptosis-alpha) (IPIA-alpha) (nm23-H6)
Other Name:	
Human Swiss Prot No.:	O75414
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of NME6 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year