

Antibodies Product Datasheet

Product Name:	BCAS2 Antibody
Catalog Number:	C14801
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BCAS2. AA range:176-225
Concentration:	1 mg/ml
Gene Name:	BCAS2
Protein Name:	Pre-mRNA-splicing factor SPF27
Other Name:	BCAS2; DAM1; Pre-mRNA-splicing factor SPF27; Breast carcinoma-amplified sequence 2; DNA amplified in mammary carcinoma 1 protein; Spliceosome-associated protein SPF 27
Human Swiss Prot No.:	O75934
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BCAS2 Polyclonal Antibody detects endogenous levels of BCAS2 protein.
Storage and Stability:	-20°C/1 year