

Antibodies Product Datasheet

Product Name:	RASL10B Antibody
Catalog Number:	C18206
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from RASL10B . at AA range: 80-160
Concentration:	1 mg/ml
Gene Name:	RASL10B
Protein Name:	Ras-like protein family member 10B
Other Name:	RASL10B; Ras-like protein family member 10B; Ras-like protein VTS58635; Ras-related protein 17; RRP17
Human Swiss Prot No.:	Q96S79
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	RASL10B Polyclonal Antibody detects endogenous levels of RASL10B protein.
Storage and Stability:	-20°C/1 year