

Antibodies Product Datasheet

Product Name:	IkappaB-alpha Antibody
Catalog Number:	R12-2938
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	IkappaB-alpha Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IkappaB-alpha.
Concentration:	1 mg/ml
Gene Name:	NFKBIA
Other Name:	I-kappa-B-alpha; IkappaBalph; IKBA; MAD3; Major histocompatibility complex enhancer-binding protein MAD3; NF-kappaB inhibitor alpha; NFKBI; NFKBIA; RL/IF-1
Human Swiss Prot No.:	P25963
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	IkappaB-alpha Antibodyantibody detects endogenous levels of IkappaB-alpha.
Storage and Stability:	Stable at -20°C for at least 1 year.