

Antibodies Product Datasheet

Product Name:	CD28 Antibody
Catalog Number:	C30340
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from T-cell-specific surface glycoprotein CD28 at AA range: 11-60
Concentration:	1 mg/ml
Gene Name:	CD28
Protein Name:	T-cell-specific surface glycoprotein CD28
Other Name:	CD28; T-cell-specific surface glycoprotein CD28; TP44; CD28
Human Swiss Prot No.:	P10747
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD28 Polyclonal Antibody detects endogenous levels of CD28 protein.
Storage and Stability:	-20°C/1 year