

Antibodies Product Datasheet

Product Name:	POP7 Antibody
Catalog Number:	C49460
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human POP7
Concentration:	1 mg/ml
Gene Name:	POP7 RPP20
Protein Name:	Ribonuclease P protein subunit p20 (RNaseP protein p20) (EC 3.1.26.5) (Ribonucleases P/MRP protein subunit POP7 homolog) (hPOP7)
Other Name:	
Human Swiss Prot No.:	O75817
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of POP7 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year