
Antibodies Product Datasheet

Product Name:	CD19 Antibody
Catalog Number:	C30333
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from B-lymphocyte antigen CD19 at AA range: 191-240
Concentration:	1 mg/ml
Gene Name:	CD19
Protein Name:	B-lymphocyte antigen CD19
Other Name:	CD19; B-lymphocyte antigen CD19; B-lymphocyte surface antigen B4; Differentiation antigen CD19; T-cell surface antigen Leu-12; CD19
Human Swiss Prot No.:	P15391
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD19 Polyclonal Antibody detects endogenous levels of CD19 protein.
Storage and Stability:	-20°C/1 year