

Antibodies Product Datasheet

Product Name:	TRIB3 Antibody
Catalog Number:	C42036
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	TRIB3 C20orf97 NIPK SKIP3 TRB3
Protein Name:	Tribbles homolog 3 (TRB-3) (Neuronal cell death-inducible putative kinase) (SINK) (p65-interacting inhibitor of NF-kappa-B)
Other Name:	
Human Swiss Prot No.:	Q96RU7
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	TRIB3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year