

## Antibodies Product Datasheet

<b>Product Name:</b>	PRX III Antibody
<b>Catalog Number:</b>	R12-2323
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000.IHC-p:1:50-300 ELISA: 1/20000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human PRX III. AA range:44-93
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	PRDX3
<b>Protein Name:</b>	Thioredoxin-dependent peroxide reductase mitochondrial
<b>Other Name:</b>	PRDX3; AOP1; Thioredoxin-dependent peroxide reductase; mitochondrial; Antioxidant protein 1; AOP-1; HBC189; Peroxiredoxin III; Prx-III; Peroxiredoxin-3; Protein MER5 homolog
<b>Human Swiss Prot No.:</b>	P30048
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	PRX III Polyclonal Antibody detects endogenous levels of PRX III protein.
<b>Storage and Stability:</b>	-20°C/1 year