

Antibodies Product Datasheet

Product Name:	Rabphilin-3A (Phospho-Ser237) Antibody
Catalog Number:	A0563
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Rabphilin 3A around the phosphorylation site of Ser237. AA range:203-252
Concentration:	1 mg/ml
Gene Name:	RPH3A
Protein Name:	Rabphilin-3A
Other Name:	RPH3A; KIAA0985; Rabphilin-3A; Exophilin-1
Human Swiss Prot No.:	Q9Y2J0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Rabphilin-3A (S237) Polyclonal Antibody detects endogenous levels of Rabphilin-3A protein only when phosphorylated at S237.
Storage and Stability:	-20°C/1 year