

Antibodies Product Datasheet

Product Name:	hCAP-G Antibody
Catalog Number:	C15234
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NCAPG. AA range:951-1000
Concentration:	1 mg/ml
Gene Name:	NCAPG
Protein Name:	Condensin complex subunit 3
Other Name:	NCAPG; CAPG; NYMEL3; Condensin complex subunit 3; Chromosome-associated protein G; Condensin subunit CAP-G; hCAP-G; Melanoma antigen NY-MEL-3; Non-SMC condensin I complex subunit G; XCAP-G homolog
Human Swiss Prot No.:	Q9BPX3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	hCAP-G Polyclonal Antibody detects endogenous levels of hCAP-G protein.
Storage and Stability:	-20°C/1 year